

FORMED ideas in motion

NEW!



ARTROMOT®-K1

CPM therapy for knee and hip joint

ARTROMOT®-K1

CPM therapy for knee and hip joint

Making good things great!

ARTROMOT®-K1 combines the optimized features of the established ARTROMOT® CPM knee devices with the newest cutting-edge technology.

ARTROMOT®-K1 has been especially developed for rental business, home care and the demanding clinical environment.

	ARTROMOT®-K1
Programming	Intuitive with pictograms
Stability	Metal frame architecture, maximum stability
Range of patient's size	120 – 200 cm
Weight	11 Kg
Range of motion	-10° / 0° / 120° knee 0° / 7° / 69° hip
Anatomical correct motion	Yes

ARTROMOT®-K1 is available in different versions:

Standard: Offering the basics of knee CPM therapy.

Comfort: Offering additional special functions, including a therapy documentation and the option to combine ARTROMOT®-K1 with electrotherapy. (Ormed's ARTROSTIM®-FOCUSplus device)

Chipcard: Stores pre-programmed treatment parameters and the therapy documentation for each single patient.

Indications, e.g. :

Operatively treated fractures and pseudarthrosis
Total knee and total hip replacement • Joint mobilization under narcosis • ACL- / PCL-reconstructions • Stable osteosyntheses
Meniscectomy • Arthroscopy • Arthrotomy

